



# Synecta™ Cell-Derived Nanoparticle Platform Offers Versatile Use of Membrane-Bound and Soluble Cytokines for Rapid and Efficient CAR-T Cell Production

Deeksha Muthumani, Chiquita Hanindya, Victor Carpio, Gil Joseph, Peter Keller, and Petr Ilyinskii



## Background

CAR-T cell manufacturing is limited by slow production and inefficient T cell activation. Current methods introduce variability and complicate workflows while not fully addressing a variety of additional technical challenges. Rapid, robust, and consistent T cell activation technology is therefore essential for advancing CAR-T manufacturing.




We developed the Synecta™ Cell-Derived Nanoparticle (CDNP) platform, which entails the surface display of membrane-bound (mb) stimulatory signals and cytokines, thus enabling more consistent, robust, and simplified CAR-T manufacturing. We hypothesized that, compared to soluble cytokines, membrane-bound formats promote more effective expansion of long-lived and highly functional T cells.

### Radically simplified process


Day 0: Synecta T1 + Transduction  
Day 3 - 9: Harvest

## Synecta Activation Platform



## Methods


Synecta T1 or T2 in combination with Synecta CK 1 or soluble cytokines were evaluated for *in vitro* T cell expansion and activation, as well as their effect on T cell phenotype.



- Isolation:** T cells were negatively isolated from healthy and patient donors.
- Activation:** T cells were activated with Synecta or comparators on Day 0.
- Transduction:** T cells were transduced with CD19 CAR LVV on Day 0 for Synecta T1 and on Day 1 for comparators.
- Cytokines:** Comparator cultures were supplemented with soluble IL-7 and IL-15 (5 and 10 ng/mL; Figure 1) or IL-2 (100 IU/mL; Figure 4). Synecta T1 was supplemented with IL-7 and IL-15 (5 and 10 ng/mL) and Synecta CK 1 (Figures 2 & 3), while Synecta T2 received IL-2 (100 IU/mL).
- Expansion:** Cells were counted on Days 3, 7, and 9 using NC-200.
- Flow Cytometry:** Activation markers (CD69, CD25) were assessed on Day 2; phenotype markers (CD4, CD8, CD45RA, CCR7, CD95) on Days 3, 7, and 9.
- Functionality:** Chronic Antigen Exposure (CAE) was performed on Day 9 post-harvest.


## Results

### Synecta T1 drives rapid, functional T-cell expansion




**Figure 1. T cell Expansion with Synecta T1 vs Commercial Activators.** (a–b) Expansion of patient T cells: Synecta T1 and Comparator 1 (n = 11) and Comparator 2 (n = 10) on Day 3 (a) and Day 9 (b). **Stats:** Mean ± SD; (a) Repeated Measures One-way ANOVA with Sidak's; (b) Wilcoxon paired t-test; \*p < 0.05, \*\*p < 0.01.

### Synecta T1 promotes robust CAR-T cell expansion without inducing exhaustion




**Figure 2. T Cell Functionality: Synecta T1 vs. Commercial Activators.** (a) Transduction efficiency of T cells activated with Synecta T1 (Day 0) versus commercial activators (Day 1), n = 2, Mean ± SD. (b) Schematic of the CAE assay. (c) CAR-T expansion comparison. (d) Synecta T1-derived CAR-T cytotoxicity at E: T 1:1 (Day 20).

### Synecta CK1 boosts Synecta T1-mediated T-cell expansion while maintaining a memory phenotype



**Figure 3. T cell expansion and Memory Phenotype with Synecta T1 supplemented with soluble and Synecta CK1 (membrane-bound) cytokines: IL-7 & IL-15.** (a) Expansion of healthy donor T cells from Day 0 to Day 9 (n = 4), with statistics shown for Day 7 (b) Memory subsets over days of expansion: T<sub>N</sub> (CD45RA<sup>+</sup>CCR7<sup>+</sup>CD95<sup>+</sup>), T<sub>SCM</sub> (CD45RA<sup>+</sup>CCR7<sup>+</sup>CD95<sup>+</sup>), T<sub>CM</sub> (CD45RA<sup>+</sup>CCR7<sup>+</sup>CD95<sup>+</sup>), T<sub>EM</sub> (CD45RA<sup>+</sup>CCR7<sup>+</sup>CD95<sup>+</sup>), T<sub>EFF</sub> (CD45RA<sup>+</sup>CCR7<sup>+</sup>CD95<sup>+</sup>) (c) CD4 and CD8 populations on Day 9 across conditions **Stats:** Mean ± SD; (a) one-way repeated measures ANOVA with Sidak's multiple comparisons; \*p < 0.05, \*\*p < 0.01.

### Synecta T2 supports versatile cytokine supplementation



**Figure 4. T-cell activation, expansion, and phenotype of the Synecta T2 compared with Comparator 1 and Synecta T1** (a) Activation indicated by CD69<sup>+</sup>CD25<sup>+</sup> double-positive cells on Day 2 across conditions. (b) Expansion in healthy donors on Day 3 (n = 4). (c) Memory phenotype (as defined in Fig. 3) on Day 3. **Stats:** Mean ± SD; (a) Friedman test with Dunn's post hoc; p < 0.05 significant, ns = not significant

## Conclusions

- Synecta CDNPs mimic APCs through membrane-bound presentation of stimulatory ligands and cytokines, eliciting synapse-like activation kinetics that promote robust T cell expansion.
- Modular design enables combining a broad range of membrane signals
- Established large-scale manufacturing at 50 L
- Three off-the-shelf products for distinct use cases:
  - Synecta T1:** Optimized activation via IL-7 and IL-15/IL-15R $\alpha$  complex
  - Synecta T2:** Flexible pairing with soluble cytokines (e.g., IL-2)
  - Synecta CK1:** Proliferation booster (e.g., for allogeneic manufacturing)
- Applicable to multiple immune cell types; expanding towards stem cell applications

